

BD Rowa™ EasyLoad

Daily use




To start your BD Rowa™ EasyLoad, press the **button next to the screen**. All required programs will be started automatically.



Place all packs to be stored on the buffer belt. To avoid dirt and foreign objects, **do not tip out the box**, but place the packs on the belt by hand.




Choose between **Stock input** or **Stock return**. Do not mix multiple stock inputs or inputs and returns to avoid inventory discrepancies.



To process stock inputs, a **delivery note number** must be entered. After entering the number, the process is confirmed with **Done**. This simplifies further processing in your pharmacy IT system.



Make sure that all packs you want to input belong to **one** order and then click the **Start input** button.



As soon as the signal lamp lights up **white**, your BD Rowa™ EasyLoad has finished the inputting process.



After each input, it is necessary to check the **disposal box**. Packs that could not be stored must be processed manually.



Packs that are output at the stock maintenance output of your robot also have to be processed **manually**.



If you want to **switch off** the BD Rowa™ EasyLoad, click the **power icon** and select **EasyLoad shut down**. The machine switches off after 30 seconds.



Want to learn more?

Visit the BD Rowa™ Customer Learning Center at <https://clc.rowa.de>.

If you have any questions, please contact **BD Rowa™ Service**.

This document does not replace the operating instructions. The contents of the operating instructions must be read carefully and followed!

bd.com/rowa

BD, the BD Logo, Rowa and EasyLoad are trademarks of Becton, Dickinson and Company or its affiliates. All other trademarks are the property of their respective owners. © 2024 BD. All rights reserved. (Jan-2024 NPM-3318)

