

# BD Rowa™ Vmotion

## Vmotion screen usage

### Hardware

#### Turn on the screen

Press the **standby button** on the supplied **remote control**.

#### NUC



In case your installation uses NUCs, they must always stay switched on. It is good practice to **restart them periodically**. Scan the QR code<sup>1</sup> to see a video on how to restart the NUC.



#### Want to learn more?


Visit the BD Rowa™ Customer Learning Center at <https://clc.rowa.de>.

### Software

#### Navigating planograms

All active planograms can be opened through the pull-down menu.

- 1 **Swipe down** in the middle of the screen to open the menu.
- 2 Switch between the **active playlist** or the **pull-down menu** using the corresponding buttons.
- 3 Use the **alphabet** to search planograms by title.
- 4 Click the **arrow** to move the planograms to the bottom.
- 5 **Tap** a planogram to open it.



You can also scroll between planograms by swiping **left** or **right**. The sequence can be defined in the BD Rowa™ Vmotion Cloud.

**Tip:** Use only **one finger** to operate the screen.

#### Restarting the software

**Swipe up** in the middle of the screen and tap the **power icon** to close the application. On the desktop you will find the BD Rowa™ Vmotion icon to start the software again.

#### Links

- 1 <https://rowa.de/en/page/training-videos-videopage/#vmotion>
- 2 <https://shop.rowa.de/en/collections/beraten-verkaufen>

If you have any questions, please contact **BD Rowa™ Service**.

This document does not replace the operating instructions. The contents of the operating instructions must be read carefully and followed!

[bd.com/rowa](https://bd.com/rowa)

BD, the BD Logo, Rowa and Vmotion are trademarks of Becton, Dickinson and Company or its affiliates. All other trademarks are the property of their respective owners. © 2024 BD. All rights reserved. (Jan-2024 NPM-3316)


#### Cleaning the screen



Use a **dry cloth** in combination with a **disinfectant cleaning spray** to clean the touchscreen without leaving streaks. Visit the BD Rowa™ Webshop<sup>2</sup> to see our range of cleaning products.

#### Product dialog

Tapping a product will open the product dialog:



#### A Product name

#### B Indication note

If you see this icon, make sure to read the notes.

#### C Alternatives and recommendations

#### D Pack sizes, prices, and stock

- ▲ : stock > 20
- ▼ : stock < 20
- ● : no stock

#### E Information, image options and videos

#### F Output product

**Tip:** Place the smallest pack size on the shelf to easily compare savings.